

Imagine... Safer, more reliable  
tissue processing



 **Histo-Tek<sup>®</sup> VP1<sup>™</sup>**  
SAKURA

## PRODUCT SPECIFICATIONS

<b>Capacity</b>	Max. 300 cassettes
<b>Dimensions (W x D x H)</b>	552 x 604 x 1.328 mm
<b>Weight</b>	130 kg
<b>Reagent reservoirs</b>	12 (including waste and cleaning)
<b>Paraffin reservoirs</b>	4
<b>Fume control</b>	Activated carbon filter and condenser
<b>Mixing</b>	Pump-in and pump-out
<b>Screen</b>	10.4 inch color LCD, touch screen
<b>Program memory</b>	10 user-defined programs
<b>Password protection</b>	3 levels



## ITEM DESCRIPTIONS

<b>1700</b>	Histo-Tek® VP1™; 100V, 50/60Hz
<b>1550</b>	Scraper, large
<b>1705</b>	Complete reagent bottle, 12 each
<b>1706</b>	Basket with handle (150 cassette size), 2 each
<b>7109</b>	Basket transport tray, 1 each
<b>6160</b>	Activated carbon filter, 2/case, 1 set

### Optional accessories

<b>6506</b>	Duct hose set (38mm OD), 1 each
<b>6507</b>	Duct hose set (75mm OD), 1 each

CONTINUOUS INNOVATION FOR PATHOLOGY



**Sakura Finetek Europe B.V.**  
The Netherlands  
Phone: +31 (0)88 592 00 00  
Fax: +31 (0)88 592 00 00 01  
sakura@sakura.eu  
www.sakura.eu